

P2RY4 Mouse mAb

Catalog No: #64043

Package Size: #64043-1 50ul #64043-2 100ul


Description

Product Name	P2RY4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human P2RY4 (AA: extra mix) expressed in E. Coli.
Target Name	P2RY4
Other Names	NRU; P2P; UNR; P2Y4
Accession No.	P51582
Calculated MW	41kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using P2RY4 mouse mAb against SW620 (1) and A549 (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using P2RY4 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using P2RY4 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.