

p63a Mouse mAb

Catalog No: #64048

Package Size: #64048-1 50ul #64048-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	p63a Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Mouse,Rat,Monkey
Immunogen Description	Synthesized peptide of human p63 α .
Target Name	p63a
Other Names	P51B; KET; p63 α
Accession No.	Q9H3D4
Calculated MW	51kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000


Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using p63 mouse mAb against A431 (1), Hela (2), Jurkat (3), THP-1 (4), NIH/3T3 (5), Cos7 (6) and PC-12 (7) cell lysate.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer (left) and lung cancer (right) using p63 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.