

PAK2 Mouse mAb

Catalog No: #64049

Package Size: #64049-1 50ul #64049-2 100ul


Description

Product Name	PAK2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of PAK2 expressed in E. Coli.
Target Name	PAK2
Other Names	PAK65; PAKgamma
Accession No.	Q13177
Calculated MW	61kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images



Western blot analysis using PAK2 mouse mAb against Hela (1), Jurkat (2), A549 (3), HEK293 (4) and K562 (5) cell lysate.



Confocal Immunofluorescence analysis of Hela cells using PAK2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human lung cancer (left) and gastric cancer (right) using PAK2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.