

PAK3 Mouse mAb

Catalog No: #64050

Package Size: #64050-1 50ul #64050-2 100ul


Description

Product Name	PAK3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of human PAK3 (AA: 1-100) expressed in E. Coli.
Target Name	PAK3
Other Names	ARA; bPAK; MRX30; MRX47; OPHN3; PAK-3; PAK3beta; beta-PAK
Accession No.	O75914
Calculated MW	62.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using PAK3 mouse mAb against Hela (1), SK-N-SH (2), T47D (3), COS7 (4), and HepG2 (5) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.