

PAPLN Mouse mAb

Catalog No: #64051

Package Size: #64051-1 50ul #64051-2 100ul


Description

Product Name	PAPLN Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human PAPLN (AA: 766-870) expressed in E. Coli.
Target Name	PAPLN
Other Names	Papilin
Accession No.	O95428
Calculated MW	137.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using PAPLN mouse mAb against HeLa (1), HepG2 (2), OCM-1 (3), Raji (4), Jurkat (5), NIH/3T3 (6) cell lysate.



Immunofluorescence analysis of NIH/3T3 cells using PAPLN mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.