

PAR4 Mouse mAb

Catalog No: #64052

Package Size: #64052-1 50ul #64052-2 100ul


Description

Product Name	PAR4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of PAR4(aa1-330) expressed in E. Coli.
Target Name	PAR4
Other Names	PAWR
Accession No.	Q96IZ0
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images




Brain



Breast carcinoma

Immunohistochemical analysis of paraffin-embedded human brain (left) and breast carcinoma (right), showing cytoplasmic and membrane localization using PAR4 mouse mAb with DAB staining.

Figure 2: Immunohistochemical analysis of paraffin-embedded human brain and breast carcinoma, showing cytoplasmic and membrane localization with DAB staining using PAWR antibody.



Immunohistochemical analysis of paraffin-embedded human Small Intestine tissues using PAR4 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.