

PAX5 Mouse mAb

Catalog No: #64056

Package Size: #64056-1 50ul #64056-2 100ul


Description

Product Name	PAX5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PAX5 (AA: 235-382) expressed in E. Coli.
Target Name	PAX5
Other Names	BSAP
Accession No.	Q02548
Calculated MW	45kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using PAX5 mouse mAb against Raji (1), and EVC-304 (2) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.