

PDGFR Mouse mAb

Catalog No: #64065

Package Size: #64065-1 50ul #64065-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PDGFR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Purified recombinant fragment of human PDGFR β expressed in E. Coli.
Target Name	PDGFR
Other Names	JTK12; PDGFR; CD140B; PDGFR1; PDGFRB
Accession No.	P09619
Calculated MW	190kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500 - 1/2000 ELISA: 1/10000

Images



Western blot analysis using PDGFR β mouse mAb against NIH/3T3 cell lysate (1).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.