

PDGFRB Mouse mAb

Catalog No: #64066

Package Size: #64066-1 50ul #64066-2 100ul


Description

Product Name	PDGFRB Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4556
Isotype	Mouse IgG1
Applications	ELISA, WB, IHC, FC
Species Reactivity	Human, Mouse
Immunogen Description	Purified recombinant fragment of human PDGFRB (AA: 33-133) expressed in E. Coli.
Target Name	PDGFRB
Other Names	IMF1; IBGC4; JTK12; PDGFR; CD140B; PDGFR1; PDGFR-1
Accession No.	P09619
Calculated MW	124kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC: 1/200 - 1/400 ELISA: 1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of NIH/3T3 cells using PDGFRB mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.