

PDK2 Mouse mAb

Catalog No: #64067

Package Size: #64067-1 50ul #64067-2 100ul


Description

Product Name	PDK2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PDK2 (AA: 178-404) expressed in E. Coli.
Target Name	PDK2
Other Names	PDHK2; PDKII
Accession No.	Q15119
Calculated MW	46.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of HeLa cells using PDK2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using PDK2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using PDK2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.