

PGA5 Mouse mAb

Catalog No: #64075

Package Size: #64075-1 50ul #64075-2 100ul


Description

Product Name	PGA5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PGA5 expressed in E. Coli.
Target Name	PGA5
Other Names	Pepsinogen 5; Pepsin A
Accession No.	P0DJD9
Calculated MW	42kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using PGA5 mouse mAb against HepG2 (1) and SMMC-7721 (2) cell lysate.



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using PGA5 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.