

PHC1 Mouse mAb

Catalog No: #64078

Package Size: #64078-1 50ul #64078-2 100ul


Description

Product Name	PHC1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PHC1 (AA: 758-1004) expressed in E. Coli.
Target Name	PHC1
Other Names	EDR1; HPH1; RAE28
Accession No.	P78364
Calculated MW	105.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HEK293 cells using PHC1 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.