

phospho-NLRC4(Ser-533) Mouse mAb

Catalog No: #64082

Package Size: #64082-1 50ul #64082-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	phospho-NLRC4(Ser-533) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Rat IgG2b
Applications	IF;FC
Species Reactivity	Human
Immunogen Description	Synthesized peptide of mouse phospho-NLRC4(Ser-533) (AA: 525-538) expressed in E. Coli.
Target Name	phospho-NLRC4(Ser-533)
Other Names	CLAN; IPAF; CLAN1; CLANA; CLANB; CLANC; CLAND; Card12; 9530011P19Rik
Accession No.	Q3UP24
Calculated MW	116.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Immunofluorescence analysis of HeLa cells using phospho-NLRC4(Ser-533) rat mAb. Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.



Immunofluorescence analysis of HeLa cells using phospho-NLRC4(Ser-533) rat mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using phospho-NLRC4(Ser-533) rat mAb (green) and negative control (red).



Flow cytometric analysis of HepG2 cells using phospho-NLRC4(Ser-533) rat mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.