

PIGR Mouse mAb

Catalog No: #64083

Package Size: #64083-1 50ul #64083-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PIGR Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PIGR (AA: extra 19-170) expressed in E. Coli.
Target Name	PIGR
Accession No.	P01833
Calculated MW	83.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using PIGR mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using PIGR mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.