

PIK3CA Mouse mAb

Catalog No: #64084

Package Size: #64084-1 50ul #64084-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PIK3CA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;IF;FC
Species Reactivity	Human,Mouse,Rat,Rabbit
Target Name	PIK3CA
Other Names	PI3K; p110-alpha
Accession No.	P42336
Calculated MW	124kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Flow cytometric analysis of HELA cells using PIK3CA mouse mAb (green) and negative control (red).




Flow cytometric analysis of HL-60 cells using PIK3CA mouse mAb (green) and negative control (red).




Immunohistochemical analysis of paraffin-embedded human cervical cancer showing cytoplasmic localization using PIK3CA mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded human colon cancer showing cytoplasmic localization using PIK3CA mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded mouse spleen showing cytoplasmic localization using PIK3CA mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded mouse kidney showing cytoplasmic localization using PIK3CA mouse mAb with DAB staining.




Immunohistochemical analysis of paraffin-embedded rat spleen showing cytoplasmic localization using PIK3CA mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded rat kidney showing cytoplasmic localization using PIK3CA mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded rabbit spleen showing cytoplasmic localization using PIK3CA mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded rabbit kidney showing cytoplasmic localization using PIK3CA mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.