

PLAGL1 Mouse mAb

Catalog No: #64092

Package Size: #64092-1 50ul #64092-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PLAGL1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PLAGL1 (AA: 118-222) expressed in E. Coli.
Target Name	PLAGL1
Other Names	ZAC; LOT1; ZAC1
Accession No.	Q9UM63
Calculated MW	50.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using PLAGL1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.