

PLIN3 Mouse mAb

Catalog No: #64098

Package Size: #64098-1 50ul #64098-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PLIN3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PLIN3 expressed in E. Coli.
Target Name	PLIN3
Other Names	PP17; TIP47; M6PRBP1
Accession No.	O60664
Calculated MW	47kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000 ICC:1/100 - 1/500 FC:1/200 - 1/400


Images



Western blot analysis using PLIN3 mouse mAb against THP-1 (1), HepG2 (2), K562 (3) and HeLa (4) cell lysate.



Immunofluorescence analysis of HeLa cells using PLIN3 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using PLIN3 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.