

PLK1 Mouse mAb

Catalog No: #64099

Package Size: #64099-1 50ul #64099-2 100ul


Description

Product Name	PLK1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PLK1 (AA: 331-508) expressed in E. Coli.
Target Name	PLK1
Other Names	PLK; STPK13
Accession No.	P53350
Calculated MW	68kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/50ELISA:1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of HeLa cells using PLK1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.