

POMC Mouse mAb

Catalog No: #64107

Package Size: #64107-1 50ul #64107-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	POMC Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human POMC (AA: 1-150) expressed in E. Coli.
Target Name	POMC
Other Names	LPH; MSH; NPP; POC; ACTH; CLIP
Accession No.	P01189
Calculated MW	29.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC: 1/200 - 1/1000 ELISA: 1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using POMC mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.