

PPARA Mouse mAb

Catalog No: #64110

Package Size: #64110-1 50ul #64110-2 100ul


Description

Product Name	PPARA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Monkey
Immunogen Description	Purified recombinant fragment of human PPARA (AA: 1-120) expressed in E. Coli.
Target Name	PPARA
Other Names	PPAR; NR1C1; hPPAR; PPARalpha
Accession No.	Q07869
Calculated MW	52.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using PPARA mouse mAb against HepG2 (1), HEK293 (2), COS7 (7), HepG2 (4), and Jurkat (5) cell lysate.



Flow cytometric analysis of Jurkat cells using PPARA mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.