

PRAK Mouse mAb

Catalog No: #64118

Package Size: #64118-1 50ul #64118-2 100ul


Description

Product Name	PRAK Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of PRAK expressed in E. Coli.
Target Name	PRAK
Other Names	PRAK; MAPKAPK5
Accession No.	Q8IW41
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IHC:1/200 - 1/1000ELISA:1/10000

Images



Immunohistochemical analysis of paraffin-embedded human liver carcinoma (A), esophagus carcinoma (B), normal spleen tissue(C), breast carcinoma (D), showing nuclear and cytoplasmic localization using PRAK mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.