

PRDM14 Mouse mAb

Catalog No: #64119

Package Size: #64119-1 50ul #64119-2 100ul


Description

Product Name	PRDM14 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PRDM14 (AA: 4-203) expressed in E. Coli.
Target Name	PRDM14
Other Names	PFM11
Accession No.	Q9GZV8
Calculated MW	64kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using mouse mAb against HEK293 (1), A549 (2), and HepG2 (3) cell lysate.



Flow cytometric analysis of A549 cells using PRDM14 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.