

PRF1 Mouse mAb

Catalog No: #64121

Package Size: #64121-1 50ul #64121-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PRF1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4557
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human PRF1 (AA: extra 148-372) expressed in E. Coli.
Target Name	PRF1
Other Names	P1; PFP; HPLH2
Accession No.	P14222
Calculated MW	61.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



All lanes : PRF1 Mouse mAb at 1/1k dilution Lane 1 : Mouse liver lysates
Lysates/proteins at 20 µg per lane. All lanes : Goat Anti-Mouse IgG H&L (HRP) at 1/20000 dilution Predicted band size: 61 kDa Observed band size: 61 kDa Exposure time: 4 seconds



Flow cytometric analysis of Jurkat cells using PRF1 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using PRF1 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.