

PROZ Mouse mAb

Catalog No: #64128

Package Size: #64128-1 50ul #64128-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PROZ Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of PROZ expressed in E. Coli.
Target Name	PROZ
Other Names	protein Z; PZ
Accession No.	P22891
Calculated MW	45kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using PROZ mouse mAb against human plasma (1).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.