

PSA (KLK3) Mouse mAb

Catalog No: #64129

Package Size: #64129-1 50ul #64129-2 100ul


Description

Product Name	PSA (KLK3) Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of KLK3 (aa26-251) expressed in E. Coli.
Target Name	PSA (KLK3)
Other Names	APS; PSA; hK3; KLK2A1; KLK3
Accession No.	P07288
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Flow cytometric analysis of PC-3 cells using KLK3 mouse mAb (right) and negative control (left).



Immunohistochemical analysis of paraffin-embedded human prostate carcinoma tissues, showing cytoplasmic localization using KLK3 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human prostate tissues using KLK3 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.