

PSAP Mouse mAb

Catalog No: #64131

Package Size: #64131-1 50ul #64131-2 100ul


Description

Product Name	PSAP Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human PSAP (AA: 17-216) expressed in E. Coli.
Target Name	PSAP
Other Names	GLBA; SAP1; SAP2
Accession No.	P07602
Calculated MW	58.1kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200-1/1000ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using PSAP mouse mAb against HEK293 (1), C6 (2), and HT1080 (3) cell lysate.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using PSAP mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissues using PSAP mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.