

PSMA Mouse mAb

Catalog No: #64135

Package Size: #64135-1 50ul #64135-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PSMA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PSMA (AA: extra 44-177) expressed in E. Coli.
Target Name	PSMA
Other Names	FOLH1; PSM; FGCP; FOLH; GCP2; mGCP; GCPII; NAALAD1; NAALAdase
Accession No.	Q04609
Calculated MW	84.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using PSMA mouse mAb against Hela (1), MCF-7 (2), HCT116 (3), and GC-7901 (4) cell lysate.



Flow cytometric analysis of HeLa cells using PSMA mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.