

PSMB8 Mouse mAb

Catalog No: #64136

Package Size: #64136-1 50ul #64136-2 100ul


Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | PSMB8 Mouse mAb                                                                          |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | WB;IHC;IF                                                                                |
| Species Reactivity    | Human,Rat                                                                                |
| Immunogen Description | Purified recombinant fragment of human PSMB8 expressed in E. Coli.                       |
| Target Name           | PSMB8                                                                                    |
| Other Names           | JMP; LMP7; D6S216; PSMB5i; RING10; D6S216E; MGC1491                                      |
| Accession No.         | P28062                                                                                   |
| Calculated MW         | 30kDa                                                                                    |
| Formulation           | Ascitic fluid containing 0.03% sodium azide.                                             |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000


Images




Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);




Western blot analysis using PSMB8 mouse mAb against Hela (1), MCF-7 (2), A431 (3), RAJI (4), MOTL4 (5) and PC-12 (6) cell lysate.



Immunofluorescence analysis of HeLa cells using PSMB8 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Immunohistochemical analysis of paraffin-embedded human intima cancer tissues using PSMB8 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human colon tissues using PSMB8 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.