

PTH1R Mouse mAb

Catalog No: #64140

Package Size: #64140-1 50ul #64140-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PTH1R Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PTH1R expressed in E. Coli.
Target Name	PTH1R
Other Names	PFE; PTHR; PTHR1; MGC138426; MGC138452; PTH1R
Accession No.	Q03431
Calculated MW	66kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000

Images



Immunofluorescence analysis of SK-BR-3 cells using PTH1R mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human kidney tissues using PTH1R mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.