

PTK6 Mouse mAb

Catalog No: #64141

Package Size: #64141-1 50ul #64141-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PTK6 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PTK6 expressed in E. Coli.
Target Name	PTK6
Other Names	BRK; FLJ42088
Accession No.	Q13882
Calculated MW	52kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using PTK6 mouse mAb against Hela (1), A549 (2) and MCF-7 (3) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.