

PTK7 Mouse mAb

Catalog No: #64142

Package Size: #64142-1 50ul #64142-2 100ul


Description

Product Name	PTK7 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human PTK7 expressed in E. Coli.
Target Name	PTK7
Other Names	CCK4; PTK7
Accession No.	Q13308
Calculated MW	118kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using PTK7 mouse mAb against Hela (1), A431 (2), HCT116 (3), Caco2 (4), HepG2 (5) and MCF-7 (6) cell lysate.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using PTK7 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.