

PYK2 Mouse mAb

Catalog No: #64150

Package Size: #64150-1 50ul #64150-2 100ul


Description

Product Name	PYK2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of PYK2 (aa815-997) expressed in E. Coli.
Target Name	PYK2
Other Names	PKB; PTK; CAKB; FAK2; FRNK; PTK2B
Accession No.	Q14289
Calculated MW	115.8kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using PYK2 mouse mAb against Raji (1), PMA induced THP-1 (2), Jurkat (3) and Ramos (4) cell lysate.



Immunohistochemical analysis of paraffin-embedded human Tonsil tissues using PYK2 mouse mAb.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.