

Rab10 Mouse mAb

Catalog No: #64151

Package Size: #64151-1 50ul #64151-2 100ul


Description

Product Name	Rab10 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human Rab10 expressed in E. Coli.
Target Name	Rab10
Other Names	RAB10
Accession No.	P61026
Calculated MW	22.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ICC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using Rab10 mouse mAb against HeLa (1) and NIH/3T3 (2) cell lysate.



Immunofluorescence analysis of LOVO cells using Rab10 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.