

Rab25 Mouse mAb

Catalog No: #64155

Package Size: #64155-1 50ul #64155-2 100ul


Description

Product Name	Rab25 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;IF;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of Rab25 expressed in E. Coli.
Target Name	Rab25
Other Names	CATX-8; RAB25
Accession No.	P57735
Calculated MW	24kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Western blot analysis using Rab25 mouse mAb against MCF-7 (1), T47D (2) and GC7901 (3) cell lysate.



Immunofluorescence analysis of A549 cells using RAB25 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Flow cytometric analysis of NIH/3T3 cells using Rab25 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human esophagus tissues (left) and human lung cancer (right) using Rab25 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.