

RAB27A Mouse mAb

Catalog No: #64156

Package Size: #64156-1 50ul #64156-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAB27A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human RAB27A (AA: FULL(1-221)) expressed in E. Coli.
Target Name	RAB27A
Other Names	GS2; RAM; RAB27; HsT18676
Accession No.	P51159
Calculated MW	24.9kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissues using RAB27A mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using RAB27A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.