

RAB4A Mouse mAb

Catalog No: #64158

Package Size: #64158-1 50ul #64158-2 100ul


Description

Product Name	RAB4A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse,Monkey
Immunogen Description	Purified recombinant fragment of human RAB4A expressed in E. Coli.
Target Name	RAB4A
Other Names	RAB4; HRES-1/RAB4
Accession No.	P20338
Calculated MW	23.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);




Western blot analysis using RAB4A mouse mAb against Jurkat (1), HeLa (2), A549 (3), HEK293 (4), K562 (5), NIH3T3 (6), PC-12 (7), and COS7 (8) cell lysate.



Flow cytometric analysis of MCF-7 cells using RAB4A mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using RAB4A mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using RAB4A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.