

RABEP2 Mouse mAb

Catalog No: #64161

Package Size: #64161-1 50ul #64161-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	RABEP2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human RABEP2 (AA: 10-200) expressed in HEK293-6e cells supernatant.
Target Name	RABEP2
Other Names	FRA
Accession No.	Q9H5N1
Calculated MW	63.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000

Images


O.D. ELISA Result




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using RABEP2 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using RABEP2 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using RABEP2 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.