

RAD23A Mouse mAb

Catalog No: #64164

Package Size: #64164-1 50ul #64164-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAD23A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human RAD23A (AA: 1-363) expressed in E. Coli.
Target Name	RAD23A
Other Names	HR23A; HHR23A
Accession No.	P54725
Calculated MW	40kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using RAD23A mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using RAD23A mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using RAD23A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.