

RON Mouse mAb

Catalog No: #64187

Package Size: #64187-1 50ul #64187-2 100ul


Description

Product Name	RON Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human RON (aa210-320) expressed in E. Coli.
Target Name	RON
Other Names	PTK8; CD136; CDw136; MST1R
Accession No.	Q04912
Calculated MW	152kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using RON mouse mAb against HCC827 (1), HT-29 (2), HCT-116 (3) and BxPC-3 (4) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.