

ROR1 Mouse mAb

Catalog No: #64190

Package Size: #64190-1 50ul #64190-2 100ul


Description

Product Name	ROR1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human
Immunogen Description	Recombinant extracellular fragment of human ROR1 (aa30-406) fused with hlgGfc tag, expressed in HEK293 cells
Target Name	ROR1
Other Names	ROR1
Accession No.	Q01973
Calculated MW	101kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using ROR1 mouse mAb against extracellular domain of human ROR1 (aa30-423).



Confocal Immunofluorescence analysis of HEK293 cells trasfected with extracellular ROR1 (aa30-406)-hlgGfc using ROR1 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.