

RPL18A Mouse mAb

Catalog No: #64192

Package Size: #64192-1 50ul #64192-2 100ul


Description

Product Name	RPL18A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human RPL18A (AA: 50-176) expressed in E. Coli.
Target Name	RPL18A
Other Names	L18A
Accession No.	Q02543
Calculated MW	20.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using RPL18A mouse mAb against NIH3T3 (1), HEK293 (2), HL60 (3), Jurkat (4), Raji (5), MOLT4 (6), and HeLa (7) cell lysate.



Immunohistochemical analysis of paraffin-embedded kidney tissues using RPL18A mouse mAb with DAB staining

Note: This product is for in vitro research use only and is not intended for use in humans or animals.