

SDHB Mouse mAb

Catalog No: #64215

Package Size: #64215-1 50ul #64215-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SDHB Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SDHB (AA: 29-280) expressed in E. Coli.
Target Name	SDHB
Other Names	IP; SDH; CWS2; PGL4; SDH1; SDH2; SDHIP
Accession No.	P21912
Calculated MW	31.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using SDHB mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using SDHB mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.