

SEC31A Mouse mAb

Catalog No: #64216

Package Size: #64216-1 50ul #64216-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SEC31A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SEC31A (AA: 429-571) expressed in E. Coli.
Target Name	SEC31A
Other Names	ABP125; ABP130; HSPC275; HSPC334; SEC31L1
Accession No.	O94979
Calculated MW	133kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using SEC31A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.