

SERPINA Mouse mAb

Catalog No: #64220

Package Size: #64220-1 50ul #64220-2 100ul


Description

Product Name	SERPINA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SERPINA (AA: 269-419) expressed in E. Coli.
Target Name	SERPINA
Other Names	PI; A1A; AAT; PI1; A1AT; nNIF; PRO2275; alpha1AT
Accession No.	P01009
Calculated MW	46.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Flow cytometric analysis of Hela cells using SERPINA mouse mAb (green) and negative control (red).



Flow cytometric analysis of Jurkat cells using SERPINA mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using SERPINA mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using SERPINA mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.