

SETD2 Mouse mAb

Catalog No: #64225

Package Size: #64225-1 50ul #64225-2 100ul


Description

Product Name	SETD2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SETD2 (AA: 2054-2245) expressed in E. Coli.
Target Name	SETD2
Other Names	LLS; HYPB; SET2; HIF-1; HIP-1; KMT3A; HBP231; HSPC069; p231HBP
Accession No.	Q9BYW2
Calculated MW	287.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HeLa cells using SETD2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.