

SETD7 Mouse mAb

Catalog No: #64226

Package Size: #64226-1 50ul #64226-2 100ul


Description

Product Name	SETD7 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;FC
Species Reactivity	Human,Rat,Monkey
Immunogen Description	Purified recombinant fragment of human SETD7 (AA: 107-366) expressed in E. Coli.
Target Name	SETD7
Other Names	KMT7; SET7; SET9; SET7/9
Accession No.	Q8WTS6
Calculated MW	40.7kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using SETD7 mouse mAb against MCF-7 (1), Hela (2), A549 (3), COS7 (4), Jurkat (5), and PC-12 (6) cell lysate.



Flow cytometric analysis of HeLa cells using SETD7 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.