

SETDB1 Mouse mAb

Catalog No: #64228

Package Size: #64228-1 50ul #64228-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SETDB1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IF
Species Reactivity	Human,Mouse,Monkey
Immunogen Description	Purified recombinant fragment of human SETDB1 expressed in E. Coli.
Target Name	SETDB1
Other Names	ESET; KG1T; KMT1E; KIAA0067; H3-K9-HMTase4; SETDB1
Accession No.	Q15047
Calculated MW	180kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000 ICC:1/200 - 1/1000 ELISA:1/10000

Images



Western blot analysis using SETDB1 mouse mAb against MCF-7 (1), T47D (2), HEK293 (3), JURKAT (4), NIH/3T3 (5), F9 (6), RAW246.7 (7) and Cos7 (8) cell lysate.



Immunofluorescence analysis of LOVO cells using SETDB1 mouse mAb (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.