

# SH3GL1 Mouse mAb

Catalog No: #64231

Package Size: #64231-1 50ul #64231-2 100ul


## Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | SH3GL1 Mouse mAb                                                                         |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | WB                                                                                       |
| Species Reactivity    | Human,Rat                                                                                |
| Immunogen Description | Purified recombinant fragment of human SH3GL1 (AA: 12-119) expressed in E. Coli.         |
| Target Name           | SH3GL1                                                                                   |
| Other Names           | EEN; CNSA1; SH3P8; SH3D2B                                                                |
| Accession No.         | Q99961                                                                                   |
| Calculated MW         | 41.5kDa                                                                                  |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide                                         |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


## Application Details

WB:1/500 - 1/2000ELISA:1/10000

## Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using SH3GL1 mouse mAb against HepG2 (1), A549 (2), HT-29 (3), and PC-12 (4) cell lysate.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.