

SHH Mouse mAb

Catalog No: #64232

Package Size: #64232-1 50ul #64232-2 100ul


Description

Product Name	SHH Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse,Monkey
Immunogen Description	Purified recombinant fragment of human SHH (AA: 26-161) expressed in E. Coli.
Target Name	SHH
Other Names	TPT; HHG1; HLP3; HPE3; SMMCI; TPTPS; MCOPCB5
Accession No.	Q15465
Calculated MW	49.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using SHH mouse mAb against LNCaP (1), HepG2 (2), PANC-1 (3), HeLa (4), SK-N-SH (5), F9 (6), NIH3T3 (7), and COS7 (8) cell lysate.



Flow cytometric analysis of HeLa cells using SHH mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using SHH mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.