

SIRT3 Mouse mAb

Catalog No: #64236

Package Size: #64236-1 50ul #64236-2 100ul


Description

Product Name	SIRT3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human SIRT3 (AA: 155-290) expressed in E. Coli.
Target Name	SIRT3
Other Names	SIR2L3
Accession No.	Q9NTG7
Calculated MW	43.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using SIRT3 mouse mAb against PANC-1 (1) cell lysate.



Flow cytometric analysis of HeLa cells using SIRT3 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.